Airstream.

Airstream Class II Type B2 (Total Exhaust) Biological Safety Cabinet, Model AB2-45_. Shown with optional support stand.

Airstream.

A

ESCO

Class II, Type B2 (Total Exhaust) Biological Safety Cabinets


The Safety Solution for Life Science Laboratories





Energy Efficient ECM Motor

- Powdered by latest generation ECM motor MADE IN USA, that is more efficient than legacy ECM and VFD motor
- 70% Energy savings compared to AC Motor
- Stable airflow, despite building voltage fluctuations & Filter Loading







	Biosafety Cabinet	Air Quality	Filtration	Electrical Safety
Standards Compliance	NSF/ ANSI 49, USA EN 12469 , Europe CFDA YY-0569, China*	ISO 14644.1 Class 3, Worldwide JIS B9920 Class 3, Japan JIS BS 5295, Class 3, UK US Fed Std 209E, Class 1, USA	EN-1822 (H14), Europe IEST-RP-CC001.3, USA IEST-RP-CC007, USA IEST-RP-CC034.1, USA	IEC 61010-1, Worldwide EN 61010-1, Europe UL 61010-1, USA CAN / CSA-22.2, No.61010-1

*CFDA certification is exclusive to AB2 models sold in China.







Main Features

- The best value of any Type B2 (Total Exhaust) Biological Safety Cabinet in the industry.
- Less energy consumption and heat output than competing products delivers lower total cost of ownership.
- The angled front, narrow profile front grille, raised armrest and frameless sash create an ergonomic work environment.
- Single piece stainless steel internal work zone eliminates welded joints where contaminants may accumulate.
- Dual-wall construction surrounds the work zone with negative pressure plenums for maximum safety.
- Fail-safe system ensures that in case of exhaust failure, the cabinet's main fan automatically shuts down to ensure safety to the user.
- Unique Esco Dynamic Chamber[™] plenum design delivers quiet, uniform airflow.
- Long life ULPA (per IEST-RP-CC001.3) supply filter and HEPA exhaust filter for airflow.
- Esco Sentinel[™] microprocessor supervises all cabinet functions.
- Esco ISOCIDE[™] antimicrobial coating on all painted surfaces minimizes contamination.

Touchpad data entry buttons permit control settings and access to diagnostics, default settings and hierarchical menus.

SET

Color coded indicator lamps display green for primary function (fan operation); blue for secondary function (fluorescent lights and electrical outlet); and orange for caution (UV lamp ON). Programmable automatic UV light timer simplifies operation, enhances contamination control, extends UV lamp life and saves energy.

uv

A graphical interface indicates cabinet performance.

Digital read-out with alpha-numeric display indicates all input, status and alarm functions.

All functions can be user activated through touch-pad programming access; see Operations Manual.

Class II Biohazard safety Cabinet

Operator, Product / Sample and Environment Protection Read all safety-related instruction before use Test / certify this cabinet at least annually

Designed to Meet IEC 61010-1 Safety / Protection Standards ISO 9001 Quality Certified Manufacturing Environment




Sentinel Microprocessor Control System, Programmable

MENU

When programmed ON • the start-up sequence confirms status with Air Safe and local time display.
the Personal Identification Number (PIN) access restricts unauthorized adjustments.

- an airflow alarm warns of deviations from normal velocities.
- an annow alarm warns or deviations from normal velocities.

<u>Airstream.</u>



ULPA-filtered air

- Unfiltered / potentially contaminated air
- Room air / Inflow air

Cabinet Filtration System


Side capture zones

Dynamic air barrier, inflow and forward-directed downflow air converge

- Ambient air is pulled through the frontgrille to prevent contamination of the work surface and work product. The inflow does not mix with the clean air within the cabinet work zone. Inflow air travels through a return path toward the common air plenum (blower plenum) at the top of the cabinet.
- Ambient air is taken in through a prefilter at the top of the cabinet, and passes through the downflow ULPA filter, entering the work zone as laminar flow. The uniform, non-turbulent air stream protects against cross contamination within and throughout the work area.
- Near the work surface, the downflow air stream splits with a portion moving toward the front air grille, and the remainder moving

to the rear air grille. A small portion of the ULPA filtered downflow enters the intake perforations at the side capture zones (small blue arrows). The uniform, non-turbulent air stream protects against cross contamination within and throughout the work area.


- A combination of inflow and downflow air streams forms an air barrier that prevents contaminated room air from entering the work zone, and prevents work surface emissions from escaping the work zone. The downflow combined with the inflow air enters the common air plenum.
- All air in the common plenum is HEPA-filtered and exhausted via a dedicated ducting system to the external environment.



Model AB2, Class II Type B2 (Total Exhaust), Biological Safety Cabinet Technical Specifications

- 1. Exhaust HEPA filter
- 2. Blower
- 3. Downflow ULPA filter
- 4. UV Light Retrofit Kit™ provision
- 5. IV-Bar Retrofit Kit[™] provision
- 6. Service fixture Retrofit Kit™
- provisions (2 on each side wall) 7. Internal single-piece
- stainless steel work zone
- 8. Electrical/ Electronic panel
- 9. Fluorescent lamp

- 10. Esco Sentinel microprocessor control system
- 11. Tempered glass sash window
- 12. Electrical outlet retrofit kit[™] provisions (2 no's)
- 13. Single-piece stainless steel work tray



- Esco Centrifugal Fan with External Rotor Motor (left) vs. Conventional Fan with Standard Motor (right)
- Esco Airstream cabinets use a combination of high performance scroll blowers (supply) and German made ebm-papst[®] permanently lubricated, centrifugal motor/blowers with external rotor designs (exhaust).
- Selected for energy efficiency, compact design, and flat profile, the completely integrated exhaust blower assembly optimizes motor cooling, with unified rotating parts and overall component balance for smooth, quiet, vibration-free operation.
- Weight is equally distributed to all bearings to extend bearing life, transfer heat and maximize speed control.



	Note to customer: Insert e	electrical voltage number ir	nto last model number digi	ts _ when ordering				
Model		AB2-35_	AB2-45_	AB2-55	AB2-65			
Nominal Size		0.9 meters (3')			 1.8 meters (6')			
External	Without Base Stand	1035 x 811 x 1460 mm 40.7" x 39.1" x 57.5"	1340 x 811 x 1460 mm 52.8" x 39.1" x 57.5"	1645 x 811 x 1460 mm 64.8" x 39.1" x 57.5"	1950 x 811 x 1460 mm 76.8" x 39.1" x 57.5"			
Dimension (W x D x H)	With Optional Base Stand, 711mm (28") type	1035 x 811 x 2171 mm 40.7" x 39.1" x 85.5"	1340 x 811 x 2171 mm 52.8" x 39.1" x 85.5"	1645 x 811 x 2171 mm 64.8" x 39.1" x 85.5"	1870 x 811 x 2171 mm 76.8" x 39.1" x85.5"			
Internal Work Area, Dimensions (W x D x H)		970 x 585 x 670 mm 38.2" x 23.0" x 26.4"	1270 x 585 x 670 mm 50.0" x 23.0" x 26.4"	1570 x 585 x 670 mm 61.8" x 23.0" x 26.4"	1870 x 585 x 670 mm 73.6" x 23.0" x 26.4"			
Internal Work Area, Space		0.43 m ² (4.67 sq.ft)	0.58 m ² (6.2 sq.ft)	0.73 m ² (7.8 sq.ft)	0.87 m² (9.3 sq.ft)			
Tested and Working Opening		173 mm (6.8") and 198 mm (7.8")						
Average	Inflow	0.53 m/s (105 fpm) at initial se						
Airflow Velocity	Downflow	0.33 m/	ith uniformity of better than	+/- 20%				
	Inflow	320 m³ <i>l</i> h (190 cfm)	419 m³/h (248 cfm)	518 m³/h (307 cfm)	617 m³ <i>l</i> h (366 cfm)			
	Downflow	622 m³/h (366 cfm)	815 m³/h (480 cfm)	1007 m³/h (593 cfm)	1200 m³/h (707 cfm)			
Airflow Volume	Certification Exhaust (Inflow + Downflow)	942 m³/h (556 cfm)	1234 m³/h (728 cfm)	1525 m³/h (900 cfm)	1817 m³/h (1072 cfm			
	Concurrent Balance Value Exhaust Volume at corresponding Static Pressure Note: Use this for HVAC sizing*	1056 m³/h (623 cfm)	1382 m³/h (816 cfm)	1708 m³/h (1008 cfm)	2035 m³/h (1201 cfm)			
	Minimum exhaust static pressure for clean exhaust filter**	465 Pa / 1.9 in H ₂ 0	364 Pa / 1.5 in H ₂ 0	330 Pa / 1.3 in H ₂ 0	417 Pa / 1.7 in H ₂ 0			
	Static Pressure with additional 174 Pa (0.7 in H_2O) required by NSF/ANSI 49:2008 Note: Use this for HVAC sizing*	639 Pa / 2.6 in H ₂ 0	538 Pa / 2.2 in H ₂ 0	504 Pa / 2.0 in H ₂ 0	591 Pa / 2.4 in H ₂ 0			
Downflow ULP	PA Filter Typical Efficiency	>99.999% for particle size between 0.1 to 0.3 microns						
Exhaust HEPA I	Filter Typical Efficiency	>99.99% at 0.3 microns						
Sound Emission***	NSF/ANSI 49	<59 dBA	<59 dBA	<60 dBA	<60 dBA			
	EN 12469	<56 dBA	<56 dBA	<57 dBA	<57 dBA			
-luorescent Lig	ht Intensity At Zero Ambient	>1000 Lux (>93 foot candles)	>1000 Lux (>93 foot candles)	>900 Lux (>84 foot candles)	>1000 Lux (>93 foot candles)			
Cabinet	Main Body	1.5 mm (0.06") 16 gauge	electro-galvanized steel with Is	ocide white oven-baked epox	y-polyester powder-coatin			
Construction	Work Zone	1.2 mm (0.05") 18 gauge	gauge electro-galvanized steel with Isocide white oven-baked epoxy-polyester powde					
	220-240V, AC, 50Hz, 1ø	AB2-3S1	AB2-4S1	AB2-5S1	AB2-6S1			
	Cabinet Full Load Amps (FLA)	2 A	2 A	2 A	2 A			
	Optional Outlets FLA	5 A	5 A	5 A	5 A			
	Cabinet Nominal Power	277 W	292 W	330 W	340 W			
Electrical****	Cabinet BTU	945	996	1126	1160			
	110-120V, AC, 60Hz, 1ø	AB2-3S2	AB2-4S2	AB2-5S2	AB2-6S2			
	Cabinet Full Load Amps (FLA)	3.5 A	3.5 A	3.5 A	3.5 A			
	Optional Outlets FLA	5 A	5 A	5 A	5 A			
	Cabinet Nominal Power	293 W	309 W	334 W	360 W			
	Cabinet BTU	1000	1054	1140	1228			
	220-240V, AC, 60Hz, 1ø	AB2-3S3	AB2-4S3	AB2-5S3	AB2-6S3			
	Cabinet Full Load Amps (FLA)	2 A	2 A	2 A	2 A			
	Optional Outlets FLA	5 A	5 A	5 A	5 A			
	Cabinet Nominal Power	293 W	308 W	345.8 W	356 W			
Cabinet BTU		1000	1051	1180	1215			
Net Weight****		175 kg (386 lbs)	229 kg (505 lbs)	238 kg (525 lbs)	279 kg (615 lbs)			
hipping Weig	ht, Maximum****	232 kg (511 lbs)	273 kg (602 lbs)	295 kg (650 lbs)	350 kg (772 lbs)			
Shipping Dime Maximum (W	nsions, x D x H)*****	1150 x 850 x 1760 mm 45.2" x 33.5" x 69.3"	1450 x 850 x 1760 mm 57.1" x 33.5" x 69.3"	1750 x 850 x 1760 mm 68.9" x 33.5" x 69.3"	2050 x 850 x 1760 m 80.7" x 33.5" x 69.3			
Shipping Volume, Maximum*****		1.72 m ³ (61 cu.ft.)	2.17 m ³ (77 cu.ft.)	2.62 m ³ (93 cu.ft.)	3.07 m ³ (108 cu.ft.)			

This Concurrent Balance Value (CBV) Exhaust (per Pitot Duct Traverse) and Static Pressure must be used when sizing the HVAC exhaust & supply.
This minimum exhaust static pressure for clean exhaust filter should <u>not</u> be used for exhaust fan sizing, and it is listed here for comparative purpose only.
Noise reading in open field condition / anechoic chamber.
Additional voltages may be available; contact Esco for ordering information.
Cabinet only, excludes optional stand.

<u>Airstream</u>,

Accesories for AB2 Biological Safety Cabinets							
Cabinet	Stainless Steel Side Wall	AB2-3S1 2010274	AB2-4S1 2010284	AB2-5S1 2010294	AB2-6S1 2010304		
Exhaust Ducting	Anti-blowback Valve	ABBV-10P 5170352					
	Air Tight Damper	B2-Damper 5170104					
Work Zone	UV Lamp		15A 0251	UV-30A 5170255			
	IV Bars	IV-965 5170250	IV-1265 5170604	IV-1565 5170278	IV-1865 5170244		
Electrical Outlet	Direct Mounted	EO-HD 5170036					
Service Fixtures	EU SF-Gas-50 mm	SF-1G50 5170015					
	EU SF-Vacuum-50 mm	SF-1V50 5170004					
	EU SF-Nitrogen-50 mm	SF-1N50 5170012					
	EU SF-Air-50 mm	SF-1A50 5170007					
	EU SF-Water-50 mm	SF-1W50 5170009					
	EU SF-Universal-50 mm	SF-2U50 5170019					
	Cu Piping SF-Must be Factory Installed	CU-Pipe 5170026					
Support Stands, Ships Flat	Support Stand with Caster Wheels (Height 28")	SPC-3A0 Gen 2 5130155	SPC-4A0 Gen 2 5130152	SPC-5A0 Gen 2 5130162	SPC-6A0 Gen 2 5130154		
	Support Stand with Caster Wheels (Height 34")	SPC-3B0 Gen 2 5130165	SPC-4B0 Gen 2 5130166	SPC-5B0 Gen 2 5130167	SPC-6B0 Gen 2 5130168		
	Support Stand with Leveling Feet - I (Height 28")	SAL-3A0 Gen 2 5130170	SAL-4A0 Gen 2 5130134	SAL-5A0 Gen 2 5130171	SAL-6A0 Gen 2 5130172		
	Support Stand with Leveling Feet - I (Height 34")	SAL-3B0 Gen 2 5130174	SAL-4B0 Gen 2 513015	SAL-5B0 Gen 2 5130176	SAL-6B0 Gen 2 5130177		



ABBV-10P



EO-H_



B2-DAMPER



N--





UV-_A-L



SAL-_A0 Gen2



ESCO GLOBAL NETWORK 43 LOCATIONS IN 22 COUNTRIES



Regional Distribution Centers



ART Equipment Biological Safety Cabinets CO₂ Incubators Compounding Pharmacy Equipment Containment / Pharma Products **Ductless Fume Hoods** Lab Animal Research Products Laboratory Centrifuges Laboratory Fume Hoods Laboratory Ovens and Incubators Laboratory Shakers Laminar Flow Clean Benches **PCR** Cabinets PCR Thermal Cyclers Powder Weighing Balance Enclosures **Ultra-low Temperature Freezers**

The Esco Group of Companies is a global life sciences tools provider with sales in over 100 countries. The group is active in lab equipment, pharma equipment and medical devices. Manufacturing facilities are located in Asia and Europe. R&D is conducted worldwide spanning the US, Europe and Asia. Sales, service and marketing subsidiaries are located in 12 major markets including the US, UK, Singapore, Japan, China and India. Regional distribution centers are located in the US, UK, and Singapore.

Life Science • Chemical Research • Assisted Reproductive Technology (ART) • Pharmaceutical Equipment • General Equipment



Esco Micro Pte. Ltd. • 21 Changi South Street 1 • Singapore 486 777 Tel +65 6542 0833 • Fax +65 6542 6920 • mail@escoglobal.com www.escoglobal.com

Esco Technologies, Inc. • 903 Sheehy Drive, Suite F, Horsham, PA 19044, USA Tel 215-441-9661 • Fax 484-698-7757 eti.admin@escoglobal.com • www.escolifesciences.us

Esco Global Offices: Bangladesh | Cameroon | China | Denmark | Germany | Hong Kong | India | Indonesia | Italy | Japan | Lithuania | Malaysia | Philippines | Russia | Singapore | South Africa | South Korea | Taiwan | Thailand | UAE | United Kingdom | USA | Vietnam









